


# PRAIRIE LABS, INC.



**BTE (Behind the Ear)**

- Standard One Year Warranty
- One Year Loss and Damage included\*
- Size 13 battery
- Gold plated battery contacts
- Switched Volume Control
- Helpful hint booklet
- One free premium zinc air battery
- Pouch



## Charged options: (Choose up to 2 Pots)

Frequency pot

Gain pot

F.B.R. pot (feedback reduction pot)

Threshold pot (if applicable)

Output Limiter pot

Variable Compression pot (if applicable)

## BTE Analog Amplifier Choices:

### Class A Amplifier

- Up to 45 dB H.F.A. full gain
- Up to 115 dB SSPL
- 0.6 MA battery drain
- Up to 380 hours of battery life with size 13 zinc air battery

### Push Pull Amplifier

- Up to 60 dB H.F.A. full on gain
- Up to 132 dB SSPL
- Battery drain of 1.5 MA
- Battery life: up to 150 hours with size 13 zinc air battery

### Class D Amplifier

- Up to 50 dB H.F.A. full on gain
- Up to 128 dB SSPL
- Battery drain of 0.7 MA
- Battery life: up to 320 hours with size 13 zinc air battery

### Low Drain Amplifier (Class A Only)

- Up to 35 dB H.F.A. full on gain
- Up to 110 dB SSPL
- Battery drain of 0.3 MA
- Battery life: up to 716 hours with size 13 zinc air battery

### AGC Amplifier

- Up to 45 dB H.F.A. full on gain
- Up to 118 dB SSPL
- Battery drain of 0.5 MA
- Battery life: up to 460 hours with size 13 zinc air battery

### AGC/Class D Amplifier

- Up to 50 dB H.F.A. full on gain
- Up to 125 dB SSPL
- Battery drain of 0.7 MA
- Battery life: up to 460 hours with size 13 zinc air battery

Battery life may vary depending on ambient noise, temperature etc.

\* A processing fee will apply: please call us at (800)322-8238, or visit our website at <http://www.prairielabs.com> for more information.