

# REVOLUTION

## RECEIVER-IN-THE-CANAL

### 4-Channel 12-Band Digital Amplifier

### 100% Digital Signal Processing

#### Features and available options:

- **Advanced** Adaptive Feedback Rejection for Improved Stability, Higher Gains, and Larger Venting at all Frequencies
- New Layered Noise Reduction Reduces Noise When Speech is not present
- 4-Channel WDRC Compression
- Excellent Performance in Noisy Environments
- Selectable Time Constant Modes
- Look-ahead Detection for Loud Impulse Noises
- Independent Expansion in Each Channel
- A variety of Receiver Units (RU) and Dome sizes available
- 1 year warranty including Loss & Damage
- Available in 5 different case colors: Beige, Gray, Charcoal, Cocoa Brown, & Black
- 12-Band Equalizer for Fine Tuning
- Up to four Programmable Memories with Memory Identification Tones
- Selection of Directional Patterns
- Wideband Output Compression limiting
- MPO Limiting Threshold Adjustments for each Channel
- Low Battery Warning Tone Alerts User of Low Battery Level
- Optional Auto T-coil or Manual T-coil
- Optional Omni Directional Microphone

Battery Size	10A
SSPL90 Max. (db SPL)	120
HF Average OSPL90 (db SPL)	110
Peak Gain (db)	45
HF Average Full on Gain (db)	40
Frequency Range (Hz)	200-8000
Batter Life (Hours)	90


#### Electronic Digital Processing Option:

- Programmable Directional Patterns --
- |                  |                |
|------------------|----------------|
| Omni Directional | Cardioid       |
| Super-cardioid   | Hyper-cardioid |



**REVOLUTION**  
Receiver-in-the-Canal  
Mini Open Fit

Revolution Suggested Fitting Range



MADE IN THE

**U S A**

Highest Quality at Lowest Prices